

BOTTLE TOP DISPENSER



BOTTLE TOP DISPENSER

Bottle Top Dispenser combined with routine calibration of equipment yield higher quality results. It makes dispensing liquids easier and the volume is adjusted quickly and set reliably. Low error, capable of reducing time and costs in research and experiment and is fully autoclavable.

Used in Chemical Handling, Pollution Testing, Soil Testing, Environmental Laboratories, Pharmaceutical Laboratories, Research and Development, Qualitative and Quantitative Analysis.

Also known as Laboratory Bottle-Top Dispenser.

PIP41-25 BOTTLE TOP DISPENSER

It is suitable for continuously dispensing liquid with a fixed volume in laboratories and is applicable in various fields

Solvent Resistance: general acid and alkali

Material: ABS



SPECIFICATIONS

Model	PIP41-25
Volume	5-25 ml
Test Volume	5 ml
Accuracy	±0.7
Gross Weight (kg)	6.1
Dimensions in cubic meters	0.72x0.21x0.2



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com