

# PIPETTE BULB



# PIPETTE BULB

Designed for efficient one-handed operation and it is made of natural rubber. Pipette bulb eliminates the need of dangerous mouth pipetting and handling of hazardous materials.

Used in Medical, Laboratory, Research, Chemistry, Biology, Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology.

Also known as Laboratory Pipette Bulb.

## PIC51-25 RUBBER PIPETTE BULB

Our pinch and release valve system allows for one-hand operation

Can accommodate pipettes up to 25mL

Easy to use

Made from natural rubber

Chemical resistant rubber filler with 3 valves for suction

Two valve system with a removable top valve for easy cleaning and one-handed use



## SPECIFICATIONS

Model	PIC51-25
Volume	2-25 ml
Color	Red
Unit	50



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)