

# SEROLOGICAL PIPETTE



# SEROLOGICAL PIPETTE

Manufactured for accurate and rapid fluid dispensing. It is sterile, pyrogen-free, cotton plugged, and can be used for a variety of liquid handling applications. The reservoir eliminates the need for multiple reagent containers and saves bench top space.

Used in Plasmid Purification, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

Also known as Laboratory Serological Pipette.

## PIC71 SEROLOGICAL PIPETTE

Easy to use

Medical grade original Polystyrene


RNase- free, DNase- free and Non-pyrogenic

Seven types of volume 1ml, 2ml, 5ml, 10ml, 25ml, 50ml and 100ml

Bidirectional graduation improve accuracy

Seven types of color code on tips make pipette

Sterilized by Gamma irradiation



## SPECIFICATIONS


Model	PIC71-4000C1	PIC71-3200C1	PIC71-1800C	PIC71-1200C
Volume	1 ml		2 ml	
Precision	2 %			
Negative Graduation (ml)	-0.4			
Graduation	1 / 100			
Color	Yellow		Green	
Sterile	Gamma			
Package	Individual	Bulk pack		
Pcs/bag	100 pcs/bag		50 pcs/bag	25 pcs/bag
Pcs/ctn	4000 pcs/carton	3200 pcs/carton	1800 pcs/carton	1200 pcs/carton
Unit	1			

Model	PIC71-0800C	PIC71-0600C	PIC71-4000C2	PIC71-3200C2
Volume	5 ml		10 ml	
Precision	2 %			
Negative Graduation (ml)	-2.5		-3	
Graduation	1 / 10			
Color	Blue		Red	
Sterile	Gamma			
Package	Individual	Bulk pack	Individual	Bulk pack
Pcs/bag	25 pcs/bag			
Pcs/ctn	800 pcs/carton	600 pcs/carton	4000 pcs/carton	3200 pcs/carton


Unit	1
------	---

Model	PIC71-3000C	PIC71-2100C	PIC71-1400C	PIC71-1000C
Volume	25 ml		50 ml	
Precision	2 %			
Negative Graduation (ml)	-3			
Graduation	1 / 5		1	
Color	Purple		Black	
Sterile	Gamma			
Package	Individual	Bulk pack	Individual	Bulk pack
Pcs/bag	25 pcs/bag			
Pcs/ctn	3000 pcs/carton	2100 pcs/carton	1400 pcs/carton	1000 pcs/carton
Unit	1			


Model	PIC71-0300I	PIC71-0300B
Volume	100 ml	
Precision	2 %	
Negative Graduation (ml)	-	-
Graduation	1	
Color	Black	
Sterile	Gamma	
Package	Individual	Bulk pack
Pcs/bag	1 pc/bag	
Pcs/ctn	300 pcs/carton	
Unit	1	




PIC71-4000C1




PIC71-3200C1




PIC71-1800C




PIC71-1200C




PIC71-0800C




PIC71-0600C




PIC71-4000C2




PIC71-3200C2




PIC71-3000C




PIC71-2100C




PIC71-1400C



PIC71-1000C



PIC71-0300I



PIC71-0300B



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)