

# UV STERILIZATION CABINET



# UV STERILIZATION CABINET

UV Sterilization Cabinet are used to sterile various equipments.They are manufactured to the highest standards to ensure reliable, sterile results.

Used in It is suitable for sterilizing high temperature resistant appliances, glass sugar porcelain cups, ceramics, production supplies, currency,non-woven fabric, composite paper, cleaning appliances, tableware,towels, clothes, shoes, packaging materials and other products..

Also known as UV Disinfection Cabinet.

## AP71 UV STERILIZATION CABINET



### SPECIFICATIONS

Model	AIP71-O253	AIP71-W253
Type	With Ozone	Without Ozone
UV Ray Wavelength	253.7/184.9nm	253.7nm
Sterilization Way	2 sides of UV ray+1 side of Ozone	3 sides of UV ray
Capacity	200 L	
Power	120 W	
UV Lamp Material	Quartz glass	
Shelves	Standard 2pc, can be optional	
Internal Size	518x523x778 mm	
External Size	750x663x1200 mm	
Packing Size	875x788x1265 mm	
Gross Weight	100 kg	



AIP71-O253



AIP71-W253





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)